

dexcomG6

MAKE KNOWLEDGE YOUR SUPERPOWER WITH DEXCOM G6



Dexcom G6 Receiver is part of the System (not shown);
smart devices sold separately.

“I’M A BIG KID
AT HEART
AND DEXCOM
HELPS ME LIVE
MY LIFE!”

Brandon D.
T1D, Dexcom Warrior



THE FUTURE OF CONTINUOUS GLUCOSE MONITORING (CGM) IS HERE

The Dexcom G6 provides real-time glucose readings every five minutes for people with type 1 or type 2 diabetes. The Dexcom G6 lets you know where your glucose is at and which direction it is heading, continuously and conveniently from your compatible smart device.*



Zero Fingersticks

Factory calibrated to allow for zero fingersticks¹



Customizable Alerts & Alarms

Set your range and get notified when you’re heading high or low



Slim, Water Resistant Wearable²

A slim water resistant sensor is discreet and easy to wear for up to 10 days



Indicated for Ages 2+ Years

The only CGM FDA-permitted for children ages 2 years and up



Smart Device Compatibility

See glucose data on your compatible smart devices in real time*



Follow app

Share your glucose data with up to 10 followers³

The Dexcom G6 requires a prescription.
*For a list of compatible devices, visit www.dexcom.com/compatibility.¹If your glucose alerts and readings from the Dexcom G6 do not match symptoms or expectations, use a blood glucose meter to make diabetes treatment decisions. ²The Dexcom G6 Sensor and Transmitter are water resistant and may be submerged under eight feet of water for up to 24 hours without failure when properly installed. ³Do not make treatment decisions based upon Dexcom Share and Dexcom Follow readings. Always confirm with your compatible smart device or receiver.

INTRODUCING THE DEXCOM G6 CONTINUOUS GLUCOSE MONITORING SYSTEM



“NO FINGERSTICKS.
BEST. NEWS. EVER.”

The Dexcom G6 eliminates fingersticks for calibration and diabetes treatment decisions.[†]

Anita V.
T1D, Dexcom Warrior



ACCURATE REAL-TIME GLUCOSE READINGS

The Dexcom G6 gives you the flexibility and convenience to monitor your glucose data when, where and how you want.

HOW DOES THE DEXCOM G6 WORK?



1 Applicator

- Simple, auto-applicator to insert the sensor just beneath the skin

2 Sensor

- Measures glucose levels for up to ten days
- Acetaminophen blocking—accurate glucose readings even when taking acetaminophen*

3 Transmitter

- Affixed to the top of sensor; sends data wirelessly to a display device

4 Display Device

Dexcom Receiver

- Touchscreen display
- Displays a glucose reading every 5 minutes – up to 288 a day
- Trend arrows and vivid colors to show the speed and direction your glucose is heading

Compatible Smart Device† (Not included with Dexcom G6 System)


- Whether at work, home, or school you can check your glucose values right on your smart device with the Dexcom G6 app

A CLOSER LOOK AT THE DEXCOM G6 MOBILE APP

The Dexcom G6 app displays glucose information at a glance. Customizable alerts and alarms can notify you of pending highs and lows.



The Dexcom G6 app is compatible with both iOS and Android devices.*



*Dexcom G6 readings can be used to make diabetes treatment decisions when taking up to a maximum acetaminophen dose of 1,000mg every 6 hours. Taking a higher dose may affect the Dexcom G6 readings.

† For a list of compatible devices, please visit dexcom.com/compatibility.

THE POWER OF SHARING

With built-in Dexcom Share, up to 10 people can follow your glucose levels, giving you the circle of support you need.*

By downloading the Dexcom Follow app†, followers can view your glucose data directly from their smart device, whether they live down the street or across the country.



Zola & Keary C.
T1Ds, Dexcom Warriors



*Do not make treatment decisions based upon Dexcom Share and Dexcom Follow readings. Always confirm with your compatible smart device or receiver.

†Internet connectivity required to access Dexcom Follow.

dexcom
CLARITY®

MANAGE YOUR DIABETES WITH DEXCOM CLARITY

Dexcom CLARITY is an important part of your Dexcom CGM system. Using your CGM data, CLARITY highlights glucose patterns, trends and statistics. You can share data with your clinic and monitor improvements between visits. Get to know Dexcom CLARITY after you set up your Dexcom CGM.

WHEN AT HOME



Log in at clarity.dexcom.com

Use your current Dexcom login or, if needed, create an account.

WHILE ON THE GO



Get weekly notifications with the Dexcom CLARITY app

Notifications available to Dexcom mobile users.





“MY A1c HAS NEVER BEEN LOWER... THAT’S THE BEAUTY OF ALWAYS KNOWING WHERE YOU ARE.”

-Vickie G.



“YOU HAVE GIVEN MY HUSBAND AND ME OUR LIVES BACK. WE’VE RETURNED TO WORK, ENROLLED THEM IN DAYCARE, AND NOW CAN SLEEP THROUGH THE NIGHT!”

-Rebeka J.



“MY DEXCOM KEEPS ME IN THE GAME!”

-Nathan P.

“DEXCOM SHARE IS LIFE-CHANGING AS SHE GOES OFF TO SCHOOL!”

-Beth R.



“THE DEXCOM HAS PROVEN TO BE AN INCREDIBLY USEFUL TOOL IN MY DIABETES CARE AND MANAGEMENT.”

-Molly J.



“NOT HAVING TO PRICK MY FINGERS FEELS AMAZING!”

-Lala J.

If your glucose alerts and readings from the G6 do not match symptoms or expectations, use a blood glucose meter to make diabetes treatment decisions.



“MY WHOLE LIFE IS ON MY PHONE... SO THE DEXCOM FITS RIGHT IN.”

-Harrel C.



“I REALLY DON’T UNDERSTAND HOW YOU CAN CONTROL DIABETES WITHOUT DEXCOM.”

-Gary K.



**Dexcom
WARRIORS**

DEXCOM SUPPORT

From paperwork to insurance processing, to product support and training, our Customer Service and Patient Support teams are here to assist you every step of the way.



DEXCOM CARE

Dexcom CGM training, software downloads, and tutorials

1-877-339-2664

Monday - Friday 6 AM - 7 PM PST (subject to change)

Saturday 6 AM - 2:30 PM PST (subject to change)

patientcare@dexcom.com



CUSTOMER SALES SUPPORT

Support with Dexcom orders and general customer questions

1-888-738-3646

Monday - Friday 6 AM - 5 PM PST

Saturday 8 AM - 2 PM PST

customerservice@dexcom.com



TECHNICAL SUPPORT

Product troubleshooting or replacement inquiries

1-844-607-8398

Available 24 hours a day; 7 days a week

techsupport@dexcom.com

Visit dexcom.com/G6 for more information.

NEXT STEPS

Ready to start the process? Follow the simple steps below to get started!

1 FILL OUT THE FORM

Please fill out the form included in your brochure and give it to your doctor's office staff to start the process to get your Dexcom.

2 LEARN MORE ABOUT DEXCOM G6

Visit dexcom.com/G6 to learn more about the benefits of Dexcom CGM.

3 TALK TO DEXCOM

A representative from Dexcom will call you at the number provided on your form. It may come up on your phone as a local number. Please listen to the message and stay on the line to speak with a Dexcom representative.

4 GET A FREE BENEFITS CHECK

Dexcom will conduct a free benefits check on your behalf and relay your out-of-pocket costs before we ship. Once you agree to your costs, we will arrange for shipment.

5 START YOUR DEXCOM

Once your boxes arrive, open your Dexcom G6 CGM System and follow the instructions. We also have video tutorials online at <http://www.dexcom.com/welcome>

If you have any questions, please call Dexcom CARE and speak with a Diabetes Coach at **1-877-339-2664**.

Dexcom®

BRIEF SAFETY STATEMENT

Failure to use the Dexcom G6 Continuous Glucose Monitoring System (G6) and its components according to the instructions for use provided with your device and available at <https://www.dexcom.com/safety-information> and to properly consider all indications, contraindications, warnings, precautions, and cautions in those instructions for use may result in you missing a severe hypoglycemia (low blood glucose) or hyperglycemia (high blood glucose) occurrence and/or making a treatment decision that may result in injury. If your glucose alerts and readings from the Dexcom G6 do not match symptoms or expectations or you're taking over the recommended maximum dosage amount of 1000mg of acetaminophen every 6 hours, use a blood glucose meter to make diabetes treatment decisions. Seek medical advice and attention when appropriate, including for any medical emergency.

Dexcom, Dexcom Follow, Dexcom CLARITY, Dexcom CARE and Dexcom G6 are registered trademarks of Dexcom, Inc. in the U.S., and may be registered in other countries.

©2019 Dexcom, Inc. All rights reserved.

LBL016046 Rev003 MT25096